



Presence-absence (P-A) broth, dehydrated 500 g

Kat. číslo produktu:

2280934

Nebezpečné

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.



Used to detect pollution in treated or potable water; described in the 16th ed., Standard Methods.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Media Type:

Presence-Absence (P-A) Broth

Veľkosť balenia:

500 g