



## Test kit, iron color disc, model IR-21

**Kat. číslo produktu:**

**2299300**

**Nebezpečné**

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Iron color disc/TPTZ low range test kit contains 100 tests.

**Simple, fast, color disc method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Technické údaje

|                            |                                  |
|----------------------------|----------------------------------|
| Model:                     | IR-21                            |
| Model testovacej súpravy:  | IR-21                            |
| Názov metódy:              | Color disc/TPTZ                  |
| Parameter:                 | Iron TPTZ                        |
| Platforma:                 | Farebný disk                     |
| Počet testov:              | 100                              |
| Rozsah merania:            | Dual Range                       |
| Rozsah merania (2):        | 0 - 1,2 mg/L Fe                  |
| Smallest increment:        | 0.01                             |
|                            | 0.05                             |
| Štandardný rozsah merania: | 0 - 2.0 mg/L Fe; 0 - 1.2 mg/L Fe |
| Štýl kufrika:              | D                                |
| Typ:                       | S dvojitou stupnicou             |