



Test kit, ascorbic acid, model ASC-1

Kat. číslo produktu:

2308100

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	ASC-1
Názov metódy:	Drop count titration/ Iodometric
Parameter:	Ascorbic Acid
Platforma:	Drop Count Titration
Počet testov:	100
Poznámka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	10 - 200 mg/L
Smallest increment:	10
Štýl kufrika:	D