



Be Right™



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Kat. číslo produktu:	2324663
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Technické údaje

Body Material:	Polycarbonate
Kapacita:	1000 mL
Veľkosť balenia:	4 /pk