



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Kat. číslo produktu:

2349500

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model:	OS-1
Model testovacej súpravy:	OS-1
Názov metódy:	Color disc/Iron reduction
Parameter:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platforma:	Farebný disk
Počet testov:	50
Prepravná hmotnosť:	4,5
Rozsah merania:	0 - 300 µg/L DEHA (a ostatní vazače kyslíku)
Rozsah merania (2):	0 - 1360 µg/L DEHA (a ostatní vazače kyslíku)
Štýl kufríka:	R