



Rack, Test Tube, Polypropylene

Kat. číslo produktu:	2407000
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Test tube rack made of Polypropylene. Will not float in water and autoclavable. Molded in one piece, stacks easily. Tube size 20 mm (0.86"); array 4 x 10. Dimensions 25 x 10.2 x 8.3 cm (9.75 x 4 x 3.25"). Nalge.

Technické údaje

Špeciálne funkcie:	Rack, Test Tube
Veľkosť balenia:	1