



Be Right™



Flask, Erlenmeyer, wide mouth, 250 mL, Subdivisions 50 mL

Kat. číslo produktu:

2489446

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included.

Technické údaje

Jednotky:	1 each
Kapacita:	250 mL
Materiál:	Glass Wide Mouth