



## Test kit, corrosion management, model CM-1

Kat. číslo produktu: 2495700

**ZASTARANÁ POLOŽKA**  
Táto položka už nie je dostupná.

To reduce lead and copper levels in distribution, service and domestic water lines, suppliers must stabilize aggressive water by adjusting pH, temperature, and hardness. Analysts can verify the success of their stabilization by screening samples with this rugged, portable test kit.

**Uses drop count titration and color disc methods**

**Maintain corrosion control of aggressive distribution system waters to meet Lead and Copper Rule Requirements**

**Establish baseline conditions to monitor treatment or disinfection strategy changes**

**Calculate Langlier and Aggressive Indexes**

---

### Technické údaje

Aplikácie:	Corrosion Management
Platforma:	Farebný disk
Počet testov:	100