



Sampler, water core, "sludge judge" large capacity, 1-1/4" diameter

Kat. číslo produktu:	2506102
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Accurate readings can be obtained on suspended solids up to a depth of 4.5 m (15'). Comes in three 1.5 m (5') sections with screw-type fittings. Designed for sampling solids in wastewater and food processing applications. Standard Sludge Judge features 3/4" diameter plastic tubing, while the large capacity Sludge Judge features 1-1/4" diameter plastic tubing for double the capacity.

Technické údaje

Priemer:	31.75 mm
Veľkosť balenia:	1 each