



Syringe filter, non-sterile PVDF, 33 mm, 0.45 µm, 50/pk

Kat. číslo produktu:	2513603
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. PVDF (polyvinylidene fluoride) filters are ideal for solutions containing proteins. Standard Female Luer-Lok inlet with standard male slip Luer outlet. Hydrophilic. Diameter: 33mm. Pore Size: 0.45µm. 50/PK.

Technické údaje

Jednotky:	50/pk
Materiál:	PVDF
Priemer:	33 mm
Veľkosť balenia:	50 /pk
Veľkosť póru:	0.45 µm