



## Parafilm M Laboratory Film

Kat. číslo produktu:	251764
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Seals labware quickly, effectively. Film prevents sample evaporation and contamination; protects anhydrous materials. Can be folded repeatedly and twisted around sharp surfaces without rupturing. Stretched lengthwise at 21 °C (71 °F), Parafilm M increases approximately 200%. Becomes more flexible, softer, stronger as temperature increases. Resistant to alcohol, medium concentration acid, and alkalies. Width: 10 cm (4"); 38 m (125') roll.

---

### Technické údaje

Veľkosť balenia:	1 each
------------------	--------