



Be Right™



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Kat. číslo produkta:

2518050

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Technické údaje

Model:	NA
Model testovacej súpravy:	AccuVac® Kit
Názov metódy:	Color disc/
	Indigo trisulfonate
Parameter:	Ozone
Platforma:	Accuvac
Počet testov:	25
Rozsah merania:	0 - 1,50 mg/L O ₃
Smallest increment:	0.05
Štandardný rozsah merania:	0 - 1.50 mg/L O ₃
Štýl kufríka:	G