



## Filter, glass fiber, diameter 47 mm, preweighed standard methods TSS procedure, 100/pk

**Kat. číslo produktu:** 2546100  
**EUR Cena bez DPH:** Kontaktujte nás  
**Dodávka do 2 týždňov**

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

---

### Technické údaje

|                  |                  |
|------------------|------------------|
| Jednotky:        | 100/pk           |
| Materiál:        | Glass Microfiber |
| Priemer:         | 47 mm            |
| Veľkosť balenia: | 100 /pk          |
| Veľkosť póru:    | 1.5 µm           |