



## Bottle rack, plastic-coated, 20 compartments

Kat. číslo produktu:	2551801
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

---

### Technické údaje

Špeciálne funkcie:	Rack, 20 Bottle Compartments
Veľkosť balenia:	1