



## Sample container, 240 mL

|                           |                 |
|---------------------------|-----------------|
| Kat. číslo produktu:      | 2552411         |
| EUR Cena bez DPH:         | Kontaktujte nás |
| Dátum odoslania neuvedené |                 |

Plastic, leak-proof. Corning Snap-Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with unique sealing mechanism so sample won't leak.

---

### Technické údaje

|                    |                |
|--------------------|----------------|
| Kapacita:          | 240 mL         |
| Priemer, vonkajší: | 80 mm (approx) |
| Veľkosť balenia:   | 140 /pk        |