



Sample transport kit

Kat. číslo produktu: 2568700

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Contains everything needed for collecting field samples, keeping them cold during transport. Includes Rubbermaid Cooler, plastic-coated rack with 15 compartments for sampling bags, 100 Whirl-Pak sampling bags (177 mL) with dechlorinating agent, one leak-proof refrigerant pack. Dimensions: 42.4 x 28.4 x 58.6 cm (16.7 x 11.2 x 23.1")