



Be Right™



pH Paper, 0 - 14 pH Range, 100 tests

Kat. číslo produktu: 2601300
EUR Cena bez DPH: Kontaktujte nás
Dostupné

Provides simple, fast, economical method of pH measurement.

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Technické údaje

Číslo metódy:	10056
Jednotky:	100 tests
Najmenší krok:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Parameter:	pH
Platforma:	Testovacie prúžky
Počet testov:	100
Rozsah merania:	0 - 14 pH units