



Cylinder, graduated, 1000 mL, w/ handles

Kat. číslo produktu:	2612953
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered” double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Technické údaje

Certifikácia:	None
Kapacita:	1000 mL
Materiál skrine:	Polypropylene
Metóda merania:	To Deliver/To Contain
Typ:	S dvojitou stupnicou
Veľkosť balenia:	1 each