



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Kat. číslo produktu:	2634001
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (Please Note: These sample vials are not autoclaveable)

Technické údaje

Dĺžka:	93mm
Jednotka:	200 /pk
Kapacita:	fill line at 100 mL (+/-2.5 mL)
Priemer:	43mm
Špeciálne funkcie:	Sample Vial, with Thiosulfate
Veľkosť balenia:	200 /pk