



## Cylinder, Graduated, 25 mL, Certified

|                           |                 |
|---------------------------|-----------------|
| Kat. číslo produktu:      | 2636240         |
| EUR Cena bez DPH:         | Kontaktujte nás |
| Dátum odoslania neuvedené |                 |

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes  $\geq 25$  mL.

---

### Technické údaje

|                  |                    |
|------------------|--------------------|
| Certifikácia:    | Class A            |
| Kapacita:        | 25 mL              |
| Materiál skrine: | Glass              |
| Metóda merania:  | To Deliver         |
| Typ:             | S jednou stupnicou |
| Veľkosť balenia: | 1 each             |
| Výška:           | 192 mm             |