



Buret, Class A, Glass, 50 mL

Kat. číslo produktu:	2636541
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Technické údaje

Kapacita:	50 mL
Materiál skrine:	Class A Glass
Špeciálne funkcie:	Buret