



Be Right™



Flask, volumetric, class A, certified, 1000 mL

Kat. číslo produktu: 2636653

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

Technické údaje

Certifikácia:	Yes
Jednotky:	1 each
Kapacita:	1000 mL
Materiál skrine:	Class A Glass