



Test kit, hypochlorite, model CN-HR

Kat. číslo produktu:

2687200

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Technické údaje

Model testovacej súpravy:	CN-HR
Názov metódy:	Drop count titration/ Thiosulfate
Parameter:	Hypochlorite, high range - (Bleach) As Cl_2
Platforma:	Drop Count Titration
Počet testov:	100
Rozsah merania:	50 - 150 g/L Cl_2 (5-15%)
Smallest increment:	0.5%
Štýl kufríka:	D