



Test strips total Hardness, 0-425 mg/L, 250 tests, individually wrapped

Kat. číslo produktu: 2793844

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Water hardness is caused almost entirely by calcium and magnesium ions. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With Hach Water Quality Test Strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Hach test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Strips are individually wrapped for maximum portability and stability. Contains 250 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Technické údaje

Jednotky:	250 tests (individually wrapped)
Model testovacej súpravy:	Test Strips
Najmenší krok:	Steps: 0, 25, 50, 120, 250, 425 ppm
Parameter:	Hardness, Total (as CaCO ₃)
Platforma:	Testovacie prúžky
Počet testov:	250
Poznámka:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Rozsah merania:	0 - 425 ppm
Rozsah merania (2):	0 - 25 gpg
Smallest increment:	Intervalos de: 0, 1.5, 3, 7, 15, 25 gpg