



Plastic ampules, pseudomonas broth, pk/50

Kat. číslo produktu:	2812250
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Pseudomonas broth plastic ampules, 2 mL. For detection of Pseudomonas species. pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Materiál:	Plastic
Media Type:	Pseudomonas Broth
Veľkosť balenia:	50 50/pk