



Plastic ampoules R2A, pk/50

Kat. číslo produktu:	2812350
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

R2A Broth Ampoules; Plastic

R2A Plastic Ampoules for detection of stressed heterotrophic bacteria; pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Technické údaje

Materiál:	Plastic
Media Type:	R2A
Veľkosť balenia:	50 /pk