



Be Right™



Bottle, Wash, Nalgene®, Narrow Mouth, 500 mL, Blue Cap/ Stem, 6/pk

Kat. číslo produktu:	2928202
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Technické údaje

Farba:	Číra
Jednotka:	6 pk
Kapacita:	500 mL
Materiál:	LDPE
Stem:	Blue Cap/Stem
Veľkosť balenia:	6 pk