



Disposable BOD bottles, 300 mL, 100/cs

Kat. číslo produktu:	2943100
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Disposable BOD Bottle

Case of 100 disposable BOD bottles without stoppers. Plastic bottle is EPA approved for BOD measurement.

Single bottles available under article number 2943101. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Fresh bottle, ready for use

No cleaning required

Unbreakable and recyclable

Oxygen will not diffuse in or out of sample

EPA compliant

Technické údaje

Body Material:	Plastic
Číslo metódy:	8043
Kapacita:	300 mL
Materiál:	PET
Veľkosť balenia:	100 100/pk