



Refill kit for 32313. CO2 calibration starter kit

Kat. číslo produktu: 32314

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Reagents to prepare a reference solution to calibrate CO2 analyzers. High carbonation 3 to 8 g/kg. Refill kit for 32313.

Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to +/- 0.5 %. Each reagent kit contains a certificate of analysis.

Quick and Easy to Use

Technické údaje

Parameter:	CO2
Podmienky skladovania:	Store in a dry place at room temperature
Prístroj:	Orbisphere 3658 Analyzer Portable CO2
Rozsah merania:	3 - 8 g/kg