



CO2 calibration starter kit, 0-3g/kg

Kat. číslo produktu: 32315

ZASTARANÁ POLOŽKA

Táto položka už nie je dostupná.

CO2 calibration starter kit, 0-3g/kg. Includes calibration bottle, adapter, and reagents from part number 32316. Suitable for ORBISPHERE instruments employing Thermal Conductivity sensors.

Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to +/- 0.5 %. Each reagent kit contains a certificate of analysis.

Quick and Easy to Use

Technické údaje

Parameter:	CO2
Podmienky skladovania:	Store in a dry place at room temperature
Prístroj:	Orbisphere 3658 Analyzer Portable CO2
Rozsah merania:	0 - 3 g/kg