



pH paper, 1.0 - 11.0 pH range

Kat. číslo produktu:	39133
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Technické údaje

Rozsah merania:	1.0 - 11.0 pH units
Veľkosť balenia:	5 4.6-meter (15') rolls per package