



MultiTest kit case (145 x 206 x 66 mm), blue polypropylene

Kat. číslo produktu:	4660100
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Technické údaje

Rozměry (V x Š x H):	145 mm x 206 mm x 66 mm
Štýl kufríka:	B
Výška:	5.75 in