



MultiTest kit case (295 x 330 x 576 mm), blue polypropylene

Kat. číslo produktu:	4660800
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Technické údaje

Rozmery (V x Š x H):	295 mm x 330 mm x 576 mm
Štýl kufrika:	J
Výška:	11 in