



OriFlo™ Filter apparatus

Kat. číslo produktu:	4966000
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Designed for rapid, qualitative filtration of particulate matter from aqueous samples using a 25 mm glass fiber filter or other medium with porosity greater than or equal to 1 micron.

Used in Hach Method 10014 for Ultra Low Range Chlorine in Wastewater

Apparatus included in the Ultra Low Range Chlorine Apparatus Set, PN 2595600

Not designed for use as a microbiological filtration device using membrane filters with a porosity of 0.45 micron or less

Cannot be autoclaved or sterilized using organic solvents