



## Beaker, griffin, low form, glass, 1000 mL

|                      |                 |
|----------------------|-----------------|
| Kat. číslo produktu: | 50053           |
| EUR Cena bez DPH:    | Kontaktujte nás |
| Dostupné             |                 |

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

---

### Technické údaje

|                  |                          |
|------------------|--------------------------|
| Jednotky:        | 1 each                   |
| Kapacita:        | 1000 mL                  |
| Materiál:        | Glass, Griffin, Low Form |
| Materiál skrine: | Glass                    |