



QL filters, pore size 8 - 12 μm , diameter 90 mm, circular, 100/ pk, qualitative

Kat. číslo produktu:	50655
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Technické údaje

Materiál:	Cellulose
Priemer:	90 mm
Veľkosť balenia:	100 /pk