



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 150 mm, circular, 100/pk, qualitative

|                           |                 |
|---------------------------|-----------------|
| Kat. číslo produktu:      | 50658           |
| EUR Cena bez DPH:         | Kontaktujte nás |
| Dátum odoslania neuvedené |                 |

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Technické údaje

|                  |           |
|------------------|-----------|
| Materiál:        | Cellulose |
| Priemer:         | 150 mm    |
| Veľkosť balenia: | 100 /pk   |