



Acid solution for EZ4004 Total Alkalinity Analyser, 0.2 N HCl (100-5000 mg/L Range), 2 L

Kat. číslo produktu:

APPC4004-04-EU **Nebezpečné**

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Technické údaje

Koncentrácia:	0.2 N HCl
Parameter:	Total Alkalinity
Veľkosť balenia:	2 L