



Digital titrator assembly

Kat. číslo produktu:	1690002
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Portable and robust digital Titrator for fast and precise titrations. High resolution (125µL/digit) and high accuracy ($\pm 1\%$). Easy handling. Applications for more than 40 parameters with interchangeable cartridges.

Accurate results - when and where you need them

Fast and Flexible

Versatile

Economical

Technické údaje

Obsahuje:	Digital Titrator device only
Vrátane príslušenstva:	Nie