



Sampling bag, sterile with dechlorinating agent, 100 mL, 100 pcs

Kat. číslo produktu:	2075333
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10mg sodium thiosulfate (will dechlorinate up to 10 mg Cl₂ in 100 mL sample) four 100mL fill-lines.

Technické údaje

Kapacita:	177 mL ml Fill Line
Veľkosť balenia:	100