



Be Right™



INTELLICAL Štandardná LDO elektróda, kábel 3 m

Kat. číslo produktu: LDO10103
EUR Cena bez DPH: Kontaktujte nás
Dostupné

Luminiscenčná technológia LDO pre laboratórium

Štandardná LDO elektróda pre digitálny multimeter, dĺžka kábla 3 m.

Testovanie v laboratóriu alebo v teréne: elektródy sú k dispozícii v štandardnej a robustnej terénnej verzii s rôznymi dĺžkami kábla.

Inteligentné a zameniteľné elektródy si uchovávajú históriu kalibrácie a umožňujú tak rýchlu výmenu elektród a jednoduché používanie s rôznymi prístrojmi. Nevyžaduje sa žiadna opakovaná kalibrácia!

Presné a správne výsledky pri nízkych i vysokých koncentráciách O₂

The Intellical LDO optical dissolved oxygen probe has no membranes to clean or replace, no electrolyte solution to replenish, and reduces calibration frequency. Trusted for dissolved oxygen measurement by thousands of facilities around the world.

Žiadna kalibrácia, žiadna polarizácia

Stored time and date stamp for each measurement, operator and sample ID, calibration history, parameter, and probe serial number

Žiadne výmeny elektrolytu, žiadne interferencie z usadenín alebo H₂S

No guesswork needed to provide the most reliable and accurate results with each measurement.

Bezchybné výsledky s minimálnym úsilím

Ideal for a multi-user environment with multiple HQD series Laboratory and Portable meters.

Technické údaje

Dĺžka:	191 mm
Dĺžka kábla:	3 m
Hmotnosť:	0,127 kg
Method Type:	Field: Luminescent Dissolved Oxygen
Parameter:	Rozpustený kyslík
Plniaci roztok:	-
Presnosť:	±0.2 mg/L for greater than 8 mg/L
Pressure accuracy:	±0.8%
Pressure Resolution:	1 mbar (1 hPa)
Product Kit:	Model: LDO101
	Accessories Included: None
Reference Junction:	NA
Rozlíšenie :	0.01/0.1 (5 digits max.)
Rozlíšenie teploty:	0.1 °C

Rozsah merania:	0.05 to 20.0 mg/L (ppm)
Saturation resolution:	0.1 %
Sensor Body Material:	Polycarbonate / ABS
Sensor Cable Length (m):	3
Sensor material:	Polycarbonate / ABS
Špeciálne funkcie:	Technológia LDO. Žiadne membrány.
Správnosť teploty:	±0.3 °C
Teplotný rozsah:	Kontinuálne použitie: 0 – 50 °C
Teplotný senzor:	# (ATC)
Test Requirements:	Parameter Needed: Dissolved Oxygen
	Minimum Sample Depth (mm): 10
Typ elektrody:	Štandardná
Typ elektródy:	Luminescent Dissolved Oxygen
Typ senzora:	Luminofor
Type of Measurement:	Field

Obsah balenia

Intellical LDO101 Laboratory Luminescent Dissolved Oxygen (LDO) probe, 3 m cable, protective shroud kit, DO cap with iButton (already mounted) Test certificate, and Basic User Manual.