



## ComLabo Počítačový softvér, pre sensION+ GLP, disk CD, kábel, adaptér USB

Kat. číslo produktu: LZW8999.99

**ZASTARANÁ POLOŽKA**  
Táto položka už nie je dostupná.

LabCom PC software, for bi-directional communication of sensION+ GLP instruments with a PC. Kit including CD, cable with 9-contact connector and USB adapter.

Complete management of the instrument from the PC

Programming and configuration of calibrations and measurements from the PC

Start calibrations and measurements from the PC

Automate the measurements of up to 4 instruments with Sampler

On line graph view of measured value

Results table personalization, the user selects the displayed information

Storage of calibration and measuring results in a data base where the data can be viewed or printed at any moment.

Print report of individual sample or results summary for several samples

Data export to Excel

### **Convenient PC Tool For Data Management and Instrument Control.**

Labcom easy allow the acquisition and the management of the obtained results from Sension+ Laboratory instruments.