



Silica Standard Solution, 500 µg/L SiO₂, 2.9 L

Kat. číslo produktu:	2100803
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Silica Standard Solution, 500 µg/L SiO₂ (2.9 L) used in the Series 5000 Silica Analyser.

Technické údaje

Koncentrácia:	500 µg/L
Parameter:	Silica
Podmienky skladovania:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Prístroj:	5500sc Silica
Veľkosť balenia:	2.9 L