



pHC3185-8 kombinovaná pH elektróda, teplotné čidlo

Kat. číslo produktu:	E16M329
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Robustá elektróda pre všeobecné aplikácie.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

Glass sensors made for long-life and durable use.

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.

Technické údaje

Materiál:	Body: epoxy
Typ elektródy:	pHC3185-8
Záruka:	6 mesiacov