



Hardness, iron, and pH test kit, model HA-62A

Kat. číslo produktu:

183701

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Aplikácie:	Water Quality
Model:	HA-62A
Parameter:	Hardness, Iron, and pH (6.6-8.4)
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	6.6 - 8.4 pH units