



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Kat. číslo produktu:

223002

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Stanovuje chlor (0-3,4 mg/L), tvrdosť (1-20 gpg), železo (0-4 mg/L) a pH (4-10).

Technické údaje

Aplikácie:	Water Quality
Model:	CN-39WR
Parameter:	Chlorine, Hardness, Iron, and pH
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0-3,4 mg/L