



**Be Right™**



## Test kit, manganese color disc, model MN-PAN

**Kat. číslo produkta:**

2350800

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

**Nebezpečné**



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Technické údaje

Model:	MN-PAN
Model testovacej súpravy:	MN-PAN
Názov metódy:	Color disc/PAN
Parameter:	Manganese
Platforma:	Farebný disk
Počet testov:	50
Prepravná hmotnosť:	3
Rozsah merania:	0 - 0,70 mg/L Mn
Smallest increment:	0.05
Štandardný rozsah merania:	0 - 0.70 mg/L Mn
Štýl kufrika:	D