



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Kat. číslo produktu:

143601

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dostupné



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Technické údaje

Model testovacej súpravy:

CA-23

Názov metódy:

Drop count titration/

Sodium hydroxide

Phenolphthalein

Parameter:

Carbon Dioxide - As CO₂

Platforma:

Drop Count Titration

Počet testov:

200

Poznámka:

*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Rozsah merania:

2 - 40 mg/L

Rozsah merania (2):

1.25 - 25 mg/L

Smallest increment:

1.25

2

5

Štýl kufríka:

D