



Be Right™



Reagent set PCB in soil pocket colorimeter II analysis system

Kat. číslo produktu: 2773500

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Reagent set for use with TPH in Soil Pocket Colorimeter Analysis System. Contains up to 18 tests.

Technické údaje

Číslo metódy:	10050
Parameter:	PCB (Polychlorinated Biphenyls)
Počet testov:	Až do 18
Prístroj:	All except DR/3000, DR/2000, DR/800
Rozsah merania:	1, 5, 10, 50 ppm thresholds