



Be Right™



Beaker, polypropylene, low form, 50 mL, 12/pk

Kat. číslo produktu: 108071
EUR Cena bez DPH: Kontaktujte nás
Vyžadat' dátum dodávky

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technické údaje

Jednotky:	12 /pk
Kapacita:	50 mL
Materiál:	Polypropylene, Low Form,
Materiál skrine:	Polypropylene
Priemer, vonkajší:	51.1 mm with spout
Volume / package size e Post Text:	/pk